



Triple Play VT® – Universal Ankle Orthosis INSTRUCTIONS FOR USE

1. Remove optional gel pack and place in freezer for a minimum of 1 hour.

NOTE: Each time the gel wrap thaws, refreeze for 1 hour.




2. Remove gel from freezer and re-attach gel pack inside ankle wrap.

3. Place a protective barrier between the gel pack and skin. A cloth sleeve is included in your package. Any thin cloth is effective. Place foot in wrap, making sure heel is comfortably positioned as shown. Secure straps around top of foot and around ankle. Make sure that

the side stay is inserted in the outside pocket to support the ankle.

4. Connect the hose to the AUX port on the Triple Play pump.
5. Press and hold the Power button until the lights on the unit are illuminated.
6. Press the MODE button on the Triple Play VT until the blue AUX light is illuminated. To increase or decrease the inflation pressure of the wrap, use the corresponding up or down arrows until the desired pressure (30 mmHg to 50 mmHg) appears on the display.

NOTE: Pressure cannot be altered during an inflation cycle.



For questions about this product call: **800.994.0464**

www.CompressionSolutions.us

OPTIONAL HAND BULB

Your product is provided with an optional hand bulb which can be used to inflate/deflate the wrap in the absence of the Triple Play VT unit. If you desire, you may use scissors to cut the tubing just below the insertion into the wrap. Then place the barbed connector on the bulb into the tubing connected to the wrap.

